



True Definition Scanner



Advanced
Technology
for Lasting
Impressions



As in most other areas of our lives, oral care is going “digital.” By choosing the 3M™ True Definition Scanner, your doctor is using digital technology to give you the best possible care. The 3M™ True Definition Scanner is a digital intraoral scanner that captures a 3D video “impression” or model of your teeth. This impression is used to produce appliances and clear aligners.

The 3M™ True Definition Scanner is:

Accurate: Resulting in nicer looking, better fitting appliances

Fast: Small, lightweight wand for fast, easy scanning

Comfortable: Especially for anyone with a low tolerance for traditional impressions

Types of Appliances

3M™ True Definition Scanner digital impressions can be used to make a wide variety of orthodontic appliances, including:

- Braces
- Clear aligners
- Spring aligners
- Retainers

Clarity™ ADVANCED Ceramic Brackets offer the ultimate combination of esthetics and performance to your patients.



What to Expect



Step 1

Your doctor or clinician simply moves the scanner around your mouth, capturing video of all sides of the teeth and gums.



Step 2

You can watch the 3D model of your teeth appear on the monitor in real time. When the scan is complete, you and your doctor can review, zoom in or out, and rotate the image to get a true picture of your oral anatomy, so you can ask questions and better understand what's being done.




Step 3

The orthodontist will use the 3D image of your teeth to determine the best treatment options for you and how they will get you to your final result.



Step 4

Should you choose a customized appliance, like customized lingual braces or aligners, your doctor will send this digital impression to a lab who will make your appliance.



“I have a very fast gag reflex which doesn’t do well when my Doctor needs to make regular impressions. The scanner makes my visits much easier and smoother!”

Kathy W., Patient



The Incognito™ Appliance System brackets are placed behind the teeth.

For more information about the 3M™ True Definition Scanner, ask your doctor or visit www.3M.com/TrueDef



Dental Products

2510 Conway Avenue
St. Paul, MN 55144-1000 USA
1-800-634-2249
Email: info3MESPE@mmm.com
Internet: www.3MESPE.com

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-888-363-3685

3M and ESPE are trademarks of
3M or 3M Deutschland GmbH. Used
under license in Canada.

Please recycle. Printed in U.S.A.
© 3M 2014. All rights reserved.
70-2013-0662-1